# reesa®

## Baby scale

**TSA0817**

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**CS**

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**DE**

**Owner's manual**

**EN**

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**FR**

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Thank you for purchasing TEESA appliance. Read this instruction manual and keep it for future reference. Producer does not take responsibility for damages caused by inappropriate handling and use of the product.

## SAFETY INSTRUCTIONS

**In order to reduce the injuries or damages, follow basic safety precautions applied when using any electrical device, including the following:**

1. Read this instruction manual carefully, even if you are familiar with the appliance. Keep the manual for future reference.
2. Place the device on flat, stable and firm surface.
3. The product is not waterproof. Do not place the appliance in water or any liquid.
4. Remove the battery if the product is not going to be used for a long time.
5. Appliance for indoor use only.
6. Do not attempt to repair this appliance yourself. Only authorized and qualified personnel may repair this device. Never disassemble this product.
7. Clean this device with soft, slightly damp cloth only. Do not use any chemical agents or abrasives to clean this product.

## BUTTONS

- UNIT key: switch units (kg / lb / st: lb)
- ON/OFF / TARE key: power on, power off / tare function

## OPERATION

- Open the battery door at the back of the scale and insert the battery (2x AAA battery). Pay attention to correct polarization. Close back the cover.
- Place the scale flat, stable surface.
- Press the "ON/OFF / TARE" button to switch on the scale. The scale is ready for operation when the display shows „0.00 kg“.
- Put the baby on the scale platform, the display shows the weight.
- Press the ON/OFF / TARE key to tare.
- Press and hold the "ON/OFF / TARE" button to switch off the scale if no more weighing is needed or it will auto off after 30 seconds of no operation.

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## TECHNICAL DATA

Maximum load: 20 kg  
 Minimal load: 50 g  
 Division: 10 g  
 Weight unit: kg / lb / st: lb  
 Material: ABS  
 Power supply: 2 x AAA battery (not included)  
 Weight: 1,15 kg  
 Platform dimensions: 50,6 x 23,3 cm  
 Device dimensions: 54,8 x 32 x 3,5 cm

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**English**  
**Correct Disposal of This Product**  
**(Waste Electrical & Electronic Equipment)**

(Applicable in the European Union and other European countries with separate collection systems) This marking shown on the product or its literature, indicates that it should not be disposed with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources. Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take this item for environmentally safe recycling. Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

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